



Electrosurgical/Electro Medical
Laparoscopic and Endoscopic Instruments




GLORIOUS ENTERPRISES

LAPAROSCOPIC SCISSORS, GRASPERS AND DISSECTORS

MONOPOLAR - SINGLE USE

- Single use with monopolar connection – offering the same characteristics as reusable instrument.
- Jaws can be easily rotatable with index finger.
- Ergonomically designed handles allowing long and fatigue free use.
- Choice of 2 handles.
- Available in a variety of jaws.
- Freedom from autoclaving.

ergo
LAPRO®



Ratchet



	ergoLAPRO	
	Non-Ratchet	Ratchet
Metzenbaum Scissors, straight, 33 cm, 5mm	500-004-33050	-
Metzenbaum Scissors, straight, 45 cm, 5mm	500-004-45050	-
Metzenbaum Scissors, curved, 33 cm, 5mm	500-014-33050	-
Metzenbaum Scissors, curved, 45 cm, 5mm	500-014-45050	-

*6 pcs in each box



	ergoLAPRO	
	Non-Ratchet	Ratchet
Endoclinch Grasping Forceps, 33 cm, 5mm	500-024-33050	501-024-33050
Endoclinch Grasping Forceps, 45 cm, 5mm	500-024-45050	501-024-45050

*6 pcs in each box



	ergoLAPRO	
	Non-Ratchet	Ratchet
Rat-Tooth Grasping Forceps, 33 cm, 5mm	500-034-33050	501-034-33050
Rat-Tooth Grasping Forceps, 45 cm, 5mm	500-034-45050	501-034-45050

*6 pcs in each box



	ergoLAPRO	
	Non-Ratchet	Ratchet
Allis Grasping Forceps, 33 cm, 5mm	500-044-33050	501-044-33050
Allis Grasping Forceps, 45 cm, 5mm	500-044-45050	501-044-45050

*6 pcs in each box



	ergoLAPRO	
	Non-Ratchet	Ratchet
Babcock Grasping Forceps, 33 cm, 5mm	500-054-33050	501-054-33050
Babcock Grasping Forceps, 45 cm, 5mm	500-054-45050	501-054-45050

*6 pcs in each box

	ergoLAPRO	
	Non-Ratchet	Ratchet
Maxi Grasping Forceps, 33 cm, 5mm	500-064-33050	501-064-33050
Maxi Grasping Forceps, 45 cm, 5mm	500-064-45050	501-064-45050
*6 pcs in each box		
	ergoLAPRO	
	Non-Ratchet	Ratchet
Oval Grasping Forceps, 33 cm, 5mm	500-074-33050	501-074-33050
Oval Grasping Forceps, 45 cm, 5mm	500-074-45050	501-074-45050
*6 pcs in each box		
	ergoLAPRO	
	Non-Ratchet	Ratchet
Long Fenestrated Forceps, 33 cm, 5mm	500-084-33050	501-084-33050
Long Fenestrated Forceps, 45 cm, 5mm	500-084-45050	501-084-45050
Short Fenestrated Forceps, 33 cm, 5mm	500-094-33050	501-094-33050
Short Fenestrated Forceps, 45 cm, 5mm	500-094-45050	501-094-45050
*6 pcs in each box		
	ergoLAPRO	
	Non-Ratchet	Ratchet
Dorsey Intestinal Forceps, 33 cm, 5mm	500-104-33050	501-104-33050
Dorsey Intestinal Forceps, 45 cm, 5mm	500-104-45050	501-104-45050
*6 pcs in each box		
	ergoLAPRO	
	Non-Ratchet	Ratchet
Atraumatic Fixation Forceps, serrated, 33 cm, 5mm	500-114-33050	501-114-33050
Atraumatic Fixation Forceps, serrated, 45 cm, 5mm	500-114-45050	501-114-45050
*6 pcs in each box		
	ergoLAPRO	
	Non-Ratchet	Ratchet
DeBakey Atraumatic Forceps, 33 cm, 5mm	500-124-33050	501-124-33050
DeBakey Atraumatic Forceps, 45 cm, 5mm	500-124-45050	501-124-45050
*6 pcs in each box		
	ergoLAPRO	
	Non-Ratchet	Ratchet
Standard Dissecting and Grasping Forceps, 33 cm, 5mm	500-134-33050	501-134-33050
Standard Dissecting and Grasping Forceps, 45 cm, 5mm	500-134-45050	501-134-45050
*6 pcs in each box		
	ergoLAPRO	
	Non-Ratchet	Ratchet
Maryland Dissector, 33 cm, 5mm	500-144-33050	501-144-33050
Maryland Dissector, 45 cm, 5mm	500-144-45050	501-144-45050
*6 pcs in each box		
	ergoLAPRO	
	Non-Ratchet	Ratchet
Duck Bill Dissector, 33 cm, 5mm	500-154-33050	501-154-33050
Duck Bill Dissector, 45 cm, 5mm	500-154-45050	501-154-45050
*6 pcs in each box		

- Individual instrument is packed in hard PET blister tray and sealed with strong and difficult to tear Tyvek lid – Offering maximum protection to sterile product.
- Method of Sterilization: Gamma radiation or ETO gas.
- *6 pcs in each inner box.



LAPAROSCOPIC ELECTRODES

MONOPOLAR - SINGLE USE

- Single use with monopolar connection
- Available with and without cable
- Freedom from autoclaving

Without Suction		
	Without Cable	With Cable
A	L-Hook Electrode, 36 cm	280-004-36050
	L-Hook Electrode, 45 cm	280-004-45050
B	J-Hook Electrode, 36 cm	280-014-36050
	J-Hook Electrode, 45 cm	280-014-45050
C	Spatula Electrode, 36 cm	280-024-36050
	Spatula Electrode, 45 cm	280-024-45050
D	Ball Electrode, 36 cm	280-034-36050
	Ball Electrode, 45 cm	280-034-45050
E	Needle Electrode, 36 cm	280-044-36050
	Needle Electrode, 45 cm	280-044-45050
F	Cone Electrode, 36 cm	280-054-36050
	Cone Electrode, 45 cm	280-054-45050
G	Knife Electrode, 36 cm	280-064-36050
	Knife Electrode, 45 cm	280-064-45050



A **B** **C** **D** **E** **F** **G**

*10 pcs in each box

*10 pcs in each box

With Suction		
	Without Cable	With Cable
A	L-Hook Electrode, 36 cm	290-004-36050
	L-Hook Electrode, 45 cm	290-004-45050
B	J-Hook Electrode, 36 cm	290-014-36050
	J-Hook Electrode, 45 cm	290-014-45050
C	Spatula Electrode, 36 cm	290-024-36050
	Spatula Electrode, 45 cm	290-024-45050
D	Ball Electrode, 36 cm	290-034-36050
	Ball Electrode, 45 cm	290-034-45050
E	Needle Electrode, 36 cm	290-044-36050
	Needle Electrode, 45 cm	290-044-45050
F	Cone Electrode, 36 cm	290-054-36050
	Cone Electrode, 45 cm	290-054-45050
G	Knife Electrode, 36 cm	290-064-36050
	Knife Electrode, 45 cm	290-064-45050



A

B

C

D


E

F

G

*10 pcs in each box

- Individual electrode is packed in sterile pouch.
- Method of Sterilization: Gamma radiation or ETO gas.
- *10 pcs in each inner box.



GLORIOUS
danmark

GLORIOUS ENTERPRISES

Naverland 11

DK 2600 Glostrup

DENMARK

TLF: +45 4362 2693

FAX: +45 4362 2611

E-mail: glorious@glorious-enterprises.dk

Homepage: www.glorious-enterprises.dk